





These Unique Strainers are Perfect For Fast and Efficient Filtering of Small Volume Samples (up to 1000 µl) Prior to FLOW or FACS Analysis.

- Preserve your sample volume; FLOWMI strainers maintain sample volume and desired cells when straining small sample sizes
- Avoid wasted time and costly delays; fast and efficient filtering using FLOWMI strainers decreases the potential for clogging of FLOW or FACS instruments
- Improves your work flow; FLOWMI strainers are oriented in a compact tray for direct press and fit attachment to a P1000 pipette tip
- Quality assurance; Sterile, FLOWMI strainers are housed in a compact, covered tray with a re-sealable bag to maintain sterility. Tray features a sliding cover for easy, one-handed use
- Use your favorite tips; FLOWMI strainers are compatible with a wide variety of P1000 tips including Fisherbrand Sure-One, VWR brand, Axygen, Nichiryo and Eppendorf

Because a Sample is a Terrible Thing to Waste

Smart. Effective. Innovative.









Press FLOWMI onto pipette tip




Strain sample into tube



Eject to dispose of tip and FLOWMI


FLOWMI™ Delivers Results You Can Count On



Cytograms of CD-Chex[™] Plus cells from left to right – Control (non-filtered); Sample filtered using Competitor; and Sample filtered using FLOWMI. Cytograms and "Percent of Parent" data demonstrate no loss of specific cell types using FLOWMI vs. competing filter.

*Flowmi is recommended for use with samples having a maximum concentration of 2 million cells/ml (40µm) or 4 million cells/ml (70µm).

FLOWMI™Ranks Best in Sample Protection



FLOWMI™ and three competing strainers were examined at four pipetted volumes (50, 200, 500, and 1000 µl). The Volume of Water Retained (µl) by each strainer, at each pipetted volume, is represented by the mean difference in weight before and after pipetting water through the strainer. (n=20 per pipetted volume; each trial using a new strainer).

Results demonstrate minimal loss of sample volume due to fluid retention inside FLOWMI™ Cell Strainer.

To see FLOWMI™ in action and for additional technical data please visit belart.com/flowmi



Cat. No.	Description
14-100-150	Sterile 40µm
03-421-228	Sterile 70µm
Patent Pending	50 Strainers /pack

Flowmi is recommended for use with samples having a maximum concentration of 2 million cells/ml (40µm) or 4 million cells/ml (70µm).

Available in 40µm and 70µm



Contact us today:

In the United States Order online: fishersci.com

Fax an order: 1-800-926-1166
Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca Fax an order: 1-800-463-2996 Call customer service: 1-800-234-7437



© 2021 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at fishersci.com/trademarks.