



DWK LIFE SCIENCES
DURAN WHEATON KIMBLE



STEVANATO GROUP EZ-fill® READY-TO-USE VIALS

Pre-sterilized Glass Vials Supporting Discovery through Delivery


Find it at fishersci.com and fishersci.ca

 **fisher scientific**
part of Thermo Fisher Scientific

STEVANATO GROUP EZ-fill®

READY-TO-USE VIALS

DWK Life Sciences offers the market's most established choice for pre-sterilized containers helping to ensure flexible and streamlined processability. As the industry shifts from a traditional, in-line vial processing model towards a more adaptable approach, pharmaceutical companies seek ways to improve efficiencies, reduce cost, and accommodate small batch or flexible filling production. Stevanato Group's EZ-fill® Ready-to-use (RTU) containment solutions allow CRO, CDMO, CMO, and Biopharma companies to maximize flexibility in aseptic filling using a common filling platform that can be integrated seamlessly into existing pharmaceutical manufacturing filling lines. EZ-fill® can improve drug packaging quality, reduce the Total Cost of Ownership (TCO), and accelerate development and time to market. EZ-fill® RTU Vials allow our customers to focus on their drug development and manufacturing competencies and to rely on DWK to meet their primary packaging needs.



KEY BENEFITS

1 FLEXIBILITY & PRODUCTION EFFICIENCIES

- Glass vials are pre-sterilized for direct entry into aseptic filling operations
- Nested configurations are designed for direct entry into flexible filling and combi line operations
- Tray configurations are suitable for small-batch and in-line filling operations

2 REDUCES COST OF OWNERSHIP


- Pre-sterilized vials supplied directly to the filling station are a cost- and space-saving alternative to operating a washing and depyrogenation line

3 REDUCED TIME TO MARKET

- The EZ-fill® pre-sterilized containment solution can be immediately filled and reduces the time from development to commercialization

4 INCREASED QUALITY & PATIENT SAFETY

- Nest & Tub and Trays formats eliminate glass to glass contact assuring the highest cosmetic quality and vial integrity
- Each vial meets exacting ISO 8362-1 Standard specifications to assure consistent fill experience and container closure integrity



STEVANATO GROUP EZ-fill®

READY-TO-USE VIALS

- Vials are designed for parenteral applications and manufactured from low extractable Type I clear borosilicate glass
- Vials are visually inspected to assure conformance to the highest dimensional and cosmetic standards
- Both configurations eliminate glass-to-glass contact to assure the highest cosmetic quality
- Manufactured with a European blowback finish

NEST & TUB

Cat. No.	Description	Dia X H (mm)	Finish (mm)	Tub Height (in)	Nest Qty	Nest/Case	Pieces/case
01-000-867	2R Vial (120), Sterile RTU Nest & Tub, Finish Up	16x35	13	2	120	15	1,800
01-000-868	4R Vial (120), Sterile RTU Nest & Tub, Finish Up	16x45	13	2	120	15	1,800
01-000-869	6R Vial (48), Sterile RTU Nest & Tub, Finish Up	22x40	20	2	48	15	720
01-000-870	10R Vial (48), Sterile RTU Nest & Tub, Finish Up	24x45	20	2	48	15	720
01-000-871	20R Vial (24), Sterile RTU Nest & Tub, Finish Up	30x55	20	3	24	15	360

PARTITIONED TRAY

Cat. No.	Description	Dia X H (mm)	Finish (mm)	Tray Qty	Tray/Case	Pieces/case
01-000-872	2R Vial (228), Sterile RTU Tray, Finish Down	16x35	13	228	10	2,280
01-000-873	4R Vial (228), Sterile RTU Tray, Finish Down	16x45	13	228	10	2,280
01-000-874	6R Vial (96), Sterile RTU Tray, Finish Down	22x40	20	96	10	960
01-000-875	10R Vial (96), Sterile RTU Tray, Finish Down	24x45	20	96	10	960
01-000-876	20R Vial (60), Sterile RTU Tray, Finish Down	30x55	20	60	8	480



NEST & TUB



PARTITIONED TRAY

ACCESSORIES & CLOSURES

TRAY ROTATION TOOL

- Tray rotation tool for vial reorientation (neck down to neck up)
- Improves workflow by assisting vial staging



Cat. No.	Description	Qty
01-258-305	Tray Rotation Tool	1

STOPPERS



- Straight plug and Igloo Stoppers
- RTU formats
- Small and large pack quantities
- Bromobutyl, Chlorobutyl
- Fluoropolymer and Silicone Coated

ALUMINUM SEALS



- Flip off, other options
- RTU formats
- Small and large pack quantities
- Ported bag options
- Blue & Red, other color options

RAYDYLYO® SNAP CAPS



- RTU Nested¹ & Bag formats
- Pre-assembled with choice of stopper
- Small and large pack quantities
- Ported bag options
- Blue & Red, other color options

¹) Compatible with Nest & Tub formats with 3" tub height.



Excellence in your hands

EZ-fill® is a registered trademark of Stevanato Group and/or its affiliates in one or more jurisdictions.
RayDyLyo® is a registered trademark of ARaymond. FBRO_030 03/2023

Distributed by Fisher Scientific. Contact us today:

In the United States

Order online: fishersci.com

Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca

Call customer service: 1-800-234-7437



© 2023 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at fishersci.com/trademarks.

BN235961555-NA