



DWK LIFE SCIENCES  
DURAN WHEATON KIMBLE

# WHEATON® COMPLETEPAK™



DESIGNED BY YOU.  
CRAFTED BY US.

Find it at [fishersci.com](https://fishersci.com) and [fishersci.ca](https://fishersci.ca)

 **fisher scientific**  
part of Thermo Fisher Scientific

# WHEATON® COMPLETEPAK

The WHEATON® **COMPLETEPAK** product line is a range of crimp top vials, rubber stoppers and seals that are ready-to-use and are off-the-shelf, which can be customized for any application. This reduces our customer's supply chain due to our one source solution capability. Each kit comes with a specific United States Pharmacopeia (USP) & European Pharmacopoeia (EP) certificate showing that the product meets or exceeds critical USP / EP standards that are enforced by the Federal Drug Administration (FDA) & European Directorate for the Quality of Medicines and Healthcare (EDQM). This kind of solution allows for our customers to focus on their core competency which is manufacturing, compounding and developing drug products.

## Target Applications

- Compounding pharmaceuticals
- Aseptic filling of injectable drug products
- Drug discovery
- Drug research
- Freeze drying / lyophilized drug products
- Large and small molecule drug packaging

## Key Benefits

- Sterilized, off-the-shelf, ready-to-use (RtU) products
- Variety of Kitted and stand-alone sterile components — a solution for any application
- Kitted solutions feature equal amounts of RtU components for reduced scrap rates
- Complete manufacturing & certificate of analysis traceability
- USP <788>/EP 2.9.19 Particulate Matter in Injections, USP <85>/EP 2.6.1, The Bacterial Endotoxins Test and USP <71>/EP 2.6.14 The Sterility Test, USP <660>/EP 3.2.1 Glass Containers for Pharmaceutical Use, and USP <381>/EP 3.2.9 Elastomeric Closures used in Parenteral Preparations
- CompletePAK Amber for light sensitive content packaging
- CompletePAK ModPAK Bulk configurations
- Ready-to-Use, Just-in-Time, minimum 180 day shelf life from time of shipment
- Container Closure Integrity (CCI) certification available upon request

## Materials

- Vials are designed for parenteral applications and manufactured from low extractable borosilicate glass
- All stoppers are manufactured from an ultrapure bromobutyl latex free formulation with extremely high chemical purity and low gas permeability
- Omniflex stoppers offer low levels of residual binding to maximize extraction of active pharmaceutical ingredients
- Omniflex stopper elasticity delivers a more effective barrier to optimize container closure integrity and content protection
- Igloo lyophilized stoppers are suited for biologics and content with sensitivities to silicone



## ENHANCING THE EFFICIENCY OF DRUG DISCOVERY AND DRUG COMPOUNDING



## TOP REASONS WHY YOU SHOULD CHOOSE COMPLETEPAK FOR YOUR APPLICATION:

1. Products supplied ready-to-use to improve workflow
2. Single source options reduce lead times & simplify administrative effort
3. Focuses your efforts on critical core competencies
4. Eliminates capital investments in sterile processing equipment
5. Frees up valuable clean room space
6. Validated processes & Certificate of Analysis traceability assures quality production outcomes
7. Saves valuable time and money while maintaining a high quality product and streamlining process


## COMPLETEPAK Kits

Clear Kit Cat. No.	Amber Kit Cat. No.	Qty/ Case	Description (pieces/kit) <i>All pieces packed sterilized</i>	ModPAK Cat. No. (Clear)	ModPAK Cat. No. (Amber)	Cases / ModPAK
06451220A	15112150A	1	2mL Vial (278), 13mm Ultra Pure Serum Stopper (285), Blue Seal (285)	15112132A	15112207	6
06451221A	15112151A	1	2mL Vial (278), 13mm Omni Flex Serum Stopper (285), Blue Seal (285)	15112133A	15112208	6
06451222A	15112152A	1	2mL Vial (278), 13mm OmniFlex Igloo Lyo Stopper (285), Blue Seal (285)	15112134A	15112209	6
06451223A	15112153A	1	2mL Vial (278), 13mm Ultra Pure Serum Stopper (285), Red Seal (285)	15112135A	15112210	6
06451224A	15112154A	1	2mL Vial (278), 13mm Omni Flex Serum Stopper (285), Red Seal (285)	15112136A	15112211	6
06451225A	15112155A	1	2mL Vial (278), 13mm OmniFlex Igloo Lyo Stopper (285), Red Seal (285)	15112137A	15112212	6
06451226A	15112156A	1	5mL Vial (120), 20mm Ultra Pure Serum Stopper (125), Blue Seal (125)	15112138A	15112213	6
06451227A	15112157A	1	5mL Vial (120), 20mm OmniFlex Serum Stopper (125), Blue Seal (125)	15112139A	15112214	6
06451228A	15112158A	1	5mL Vial (120), 20mm OmniFlex Igloo Lyo Stopper (125), Blue Seal (125)	15112140A	15112215	6
06451229A	15112159A	1	5mL Vial (120), 20mm Ultra Pure Serum Stopper (125), Red Seal (125)	15112141A	15112216	6
06451230A	15112160A	1	5mL Vial (120), 20mm OmniFlex Serum Stopper (125), Red Seal (125)	15112142A	15112217	6
06451231A	15112161A	1	5mL Vial (120), 20mm OmniFlex Igloo Lyo Stopper (125), Red Seal (125)	15112143A	15112218	6
06451232A	15112162A	1	10mL Vial (104), 20mm Ultra Pure Serum Stopper (125), Blue Seal (125)	15112144A	15112219	6
06451233A	15112163A	1	10mL Vial (104), 20mm OmniFlex Serum Stopper (125), Blue Seal (125)	15112145A	15112220	6
06451234A	15112164A	1	10mL Vial (104), 20mm OmniFlex Igloo Lyo Stopper (125), Blue Seal (125)	15112146A	15112221	6
06451235A	15112165A	1	10mL Vial (104), 20mm Ultra Pure Serum Stopper (125), Red Seal (125)	15112147A	15112222	6
06451236A	15112166A	1	10mL Vial (104), 20mm OmniFlex Serum Stopper (125), Red Seal (125)	15112148A	15112223	6
06451237A	15112167A	1	10mL Vial (104), 20mm OmniFlex Igloo Lyo Stopper (125), Red Seal (125)	15112149A	15112224	6

## COMPLETEPAK Sterile Components\*

Cat. No.	Description	Qty / Case
<b>Clear Vials</b>		
06451238A	2mL Clear Sterile Serum Vial	278
06451239A	5mL Clear Sterile Serum Vial	120
06451240A	10mL Clear Sterile Serum Vial	104
<b>Amber Vials</b>		
15112168A	2mL Amber Sterile Serum Vial	278
15112169A	5mL Amber Sterile Serum Vial	120
15112170A	10mL Amber Sterile Serum Vial	104

Cat. No.	Description	Qty / Case
<b>Stoppers</b>		
06451241A	13mm OmniFlex Sterile Serum Stopper	285
06451242A	13mm Ultra Pure Sterile Serum Stopper	285
06451243A	13mm OmniFlex Sterile Igloo Lyo Stopper	285
06451244A	20mm OmniFlex Sterile Serum Stopper	125
06451245A	20mm Ultra Pure Sterile Serum Stopper	125
06451246A	20mm OmniFlex Sterile Igloo Lyo Stopper	125
<b>Seals</b>		
06451247A	13mm Sterile Red Flip Seal	285
06451248A	13mm Sterile Blue Flip Seal	285
06451249A	20mm Sterile Red Flip Seal	125
06451250A	20mm Sterile Blue Flip Seal	125



Excellence in your hands

WHEATON® is a registered trademark of DWK Life Sciences LLC.  
Copyright © 2021 DWK Life Sciences LLC. FBRO\_018 07/2021

Distributed by Fisher Scientific. Contact us today:

**In the United States**

Order online: [fishersci.com](http://fishersci.com)

Fax an order: 1-800-926-1166

Call customer service: 1-800-766-7000

**In Canada**

Order online: [fishersci.ca](http://fishersci.ca)

Fax an order: 1-800-463-2996

Call customer service: 1-800-234-7437



© 2021 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at [fishersci.com/trademarks](http://fishersci.com/trademarks).

BN20171340